



Microbix Biosystems Inc.
 265 Watline Avenue, Mississauga
 Ontario Canada L4Z 1P3
 Tel 416-234-1624 Fax 416-234-1626

WFI Quality Water 0.1 micron filtered

Catalogue Number: WF-60-200
Product Lot #: 60039
Expiry Date: June 2013

Description:

Water for pharmaceutical purposes. This product is intended for further manufacturing or research use. Not for diagnostic, therapeutic use or parenteral administration.

Packaging: USP Class VI, Plastic Sterile Container

Storage: Store at 2°C to 30°C

Preparation:

Microbix' WFI Quality Water is prepared from Water For Injection that is 0.1µm filtered and packaged in sterile containers. It contains no antimicrobial agents or other added substances.

Testing : Tested according to USP Packaged Sterile Purified Water and Packaged Sterile Water for Injection.

Parameter	Method	Specification	Results
Appearance	Visual Inspection	Clear Liquid	Clear Liquid
Sterility	USP <71>	No Growth	No growth
Endotoxin	USP <85>	< 0.25 EU/ml	<0.03 EU/ml
Oxidizable Substances	Current USP	Pass	Pass
Conductivity	USP <645>	≤ 5 µS/cm	1.188 µS/cm
Particulate Matter for Large-Volume Injections	USP <788>	Meets USP Requirements	Meets USP requirements

 Quality Assurance

 Date