

Product Information Sheet

Microbix Biosystems Inc.

265 Watline Avenue, Mississauga
Ontario Canada L4Z 1P3
Tel 416-234-1624 Fax 416-234-1626



WFI Quality Water 0.1 micron filtered

Catalogue Number: WF-60-020

Product Lot #: 60040

Expiry Date: July 2013

Description:

Water for pharmaceutical purposes. This product is intended for further manufacturing or research use. Not for diagnostic, therapeutic use or parenteral administration.

Packaging: USP Class VI, Plastic Sterile Container

Storage: Store at 2°C to 30°C

Preparation:

Microbix' WFI Quality Water is prepared from Water For Injection that is 0.1µm filtered and packaged in sterile containers. It contains no antimicrobial agents or other added substances.

Testing : Tested according to EP Requirements

Parameter	Specifications	Results
Appearance	Clear Liquid	Clear Liquid
Sterility	No Growth	No Growth
Endotoxin	<0.25 EU/ml	<0.03 EU/ml
Oxidizable Substances	The solution remains faintly pink	Conforms
Acidity	Meets EP Requirements	Conforms
Chloride	The solution shows no change in appearance for at least 15 minutes	Conforms
Sulphate	The solution shows no change in appearance for at least 1hour	Conforms
Ammonium	Maximum 0.2 ppm	LT 0.2ppm, Conforms
Calcium and Magnesium	A pure blue colour is produced	Conforms
Residue on Evaporation	Maximum 0.001%	0.001%, Conforms
Nitrate	NMT 0.2ppm	LT 0.2ppm, Conforms
Aluminium	Maximum 10ppb	LT 10ppb, Conforms
Heavy Metals	NMT 0.1ppm	LT 0.1ppm, Conforms
Conductivity	NMT 5.1µS/cm	0.84 µS/cm, Conforms
Particulate Contamination	Meets EP Requirements	Conforms

Quality Assurance

Date