



**Microbix Biosystems Inc.**  
 265 Watline Avenue, Mississauga  
 Ontario Canada L4Z 1P3  
 Tel 416-234-1624 Fax 416-234-1626

**WFI Quality Water 0.1 micron filtered**

*Catalogue Number: WF-60-020*  
*Product Lot #: 60040*  
*Expiry Date: July 2013*

*Description:*

Water for pharmaceutical purposes. This product is intended for further manufacturing or research use. Not for diagnostic, therapeutic use or parenteral administration.

*Packaging:* USP Class VI, Plastic Sterile Container


*Storage:* Store at 2°C to 30°C


*Preparation:*

Microbix' WFI Quality Water is prepared from Water For Injection that is 0.1µm filtered and packaged in sterile containers. It contains no antimicrobial agents or other added substances.

*Testing :* Tested according to USP Packaged Sterile Purified Water and Packaged Sterile Water for Injection.

Parameter	Method	Specification	Results
Appearance	Visual Inspection	Clear Liquid	Clear Liquid
Sterility	USP <71>	No Growth	No growth
Endotoxin	USP <85>	< 0.25 EU/ml	<0.03 EU/ml
Oxidizable Substances	Current USP	Pass	Pass
Conductivity	USP <645>	≤ 5 µS/cm	0.84 µS/cm
Particulate Matter for Large-Volume Injections	USP <788>	Meets USP Requirements	Meets USP requirements

  
 \_\_\_\_\_  
 Quality Assurance

  
 \_\_\_\_\_  
 Date