



**Microbix Biosystems Inc.**  
265 Watline Avenue  
Mississauga, Ontario, Canada L4Z 1P3  
Tel 905-361-8910 Fax 905-361-8911

*Product Name:*

## Negative Control Sample

*Catalogue Number:* PT-99-01

*Storage:* Store this preparation at 2°C to 8°C.

*Description:* This product contains negative antigen control in storage buffer.

*Intended use:* The product is to be used as a negative control sample in commercially-available clinical antigen test kits, and in non-commercial clinical laboratory tests for the detection of viral antigens. The product is commonly used as a challenge sample in laboratory proficiency testing programs. It is recommended that product is tested according to kit Manufacturer's instructions for testing of clinical samples.

*Hazards:* We are aware of no specific biological hazards associated with this product. The preparation should contain no infectious material. Generally accepted good laboratory practices appropriate to biological reagents should be employed in its handling. The buffer solution contains 0.02% sodium azide (CAS No. 26628-22-8) as a preservative.

*Assistance:* If you have any questions regarding the production, testing or use of this antigen, please send them together with any relevant data to Microbix Technical Services, at [customer.service@microbix.com](mailto:customer.service@microbix.com) or fax 905-361-8911. Your complete satisfaction with the performance of this product is important to us.

## Quality Control Information

*Product name:* Negative Control Sample

*Lot number:* 250811

Test Kit	Target	Result
SAS Adeno Test Catalogue # 067020	Adenovirus	<b>Negative</b>
SAS Rota Test Catalogue # 069020	Rotavirus	<b>Negative</b>
Quick Vue Influenza A & B Catalogue # 20183	Influenza A & B	<b>Negative</b>
Xpect RSV Catalogue # R24601	Respiratory Syncytial Virus	<b>Negative</b>
Xpext Influenza A and B Catalogue # R24600	Influenza A & B	<b>Negative</b>
Quick Vue Chlamydia Catalogue # OB006	Chlamydia	<b>Negative</b>

*Quality Control Summary:* This lot meets specifications and is released for shipment.

*Expiry Date:* 09-January, 2012

Quality Assurance Signature

23-September, 2011