



Microbix Biosystems Inc.

265 Watline Avenue, Mississauga
 Ontario, Canada, L4Z 1P3
 Tel 416- 234-1624 Fax 416- 234-1626

USP Purified Water

Catalogue Number: **WF-70-020**
Product Lot #: **70009**
Retest Date: **15 June 2012**

Description:

Purified Water is intended for use as an ingredient of official preparations and in tests and assays. This product is not sterile. Do not use Purified Water for sterile dosage forms or preparations intended for parenteral administration.

Storage:

Store at 2°C to 30°C

Preparation:

USP Purified water is not filled in a sterile area and into sterilized containers. It is produced from high purity water that meets USP microbial quality requirements. It contains no added substance.

Testing:

Meets USP requirements.

Parameter	Specifications	Results
Oxidizable Substances	The pink colour doesn't disappear	Pass
Conductivity	NMT 5µS/cm	0.91 µS/cm
Total Organic Carbon	NMT 500 ppb	41 ppb
Microbial Count	NMT 100 cfu/ml (at time of manufacture)	< 1 cfu/ml

Quality Assurance

Date