

# MICROBIX'S MOST POPULAR MULTIPLEX QUALITY CONTROLS

Enhance test monitoring and Performance with Microbix's  
**Respiratory** External Quality Controls (**QAPs™**)



## FEATURES

- **Whole workflow** controls
- **Cross-platform compatible**
- **Copan FLOQSwab®** stable for **2 years**
- Room-temperature stable swabs at **2-30°C**.
- Samples eluted in your choice of media
- **Verified performance** on IVD platforms
- Inactivated to **meet clinical laboratory biosafety requirements**

## PRODUCT INTENDED USE

- **REDx™ products** are unassayed controls that can be used for ongoing **External Quality Control (QC)**, benchmarking **laboratory testing performance**, operator **training & competency** and end-to-end **workflow evaluation**.
- **PROCEEDx™ products** are used for **assay verification/validation, training, and R&D purposes**. Products are for **RESEARCH USE ONLY!**

## RESPIRATORY VIRAL PANEL

### SARS-CoV-2/Flu A & Flu B/RSV Positive Swab

Formulated with inactivated whole-genome SARS-CoV-2 cDNA and inactivated native Influenza A, Influenza B, and Respiratory Syncytial Virus

**MICROBIX BRAND  
CATALOGUE NUMBER**

**For RUO: VP-S-19-M1  
For IVD: RED-S-19-M1**

