

MICROBIX'S MOST POPULAR MULTIPLEX QUALITY CONTROLS

Enhance test monitoring and Performance with Microbix's
STI External Quality Controls (QAPs™)



FEATURES

- **Whole workflow** controls
- **Cross-platform compatible**
- **Copan FLOQSwab®** stable for **2 years**
- Room-temperature stable swabs at **2-30°C**.
- Samples eluted in your choice of media
- **Verified performance** on IVD platforms
- Inactivated to **meet clinical laboratory biosafety requirements**

PRODUCT INTENDED USE

- **REDx™ products** are unassayed controls that can be used for ongoing **External Quality Control (QC)**, benchmarking **laboratory testing performance**, operator **training & competency** and end-to-end **workflow evaluation**.
- **PROCEEDx™ products** are used for **assay verification/validation, training, and R&D purposes**. Products are for **RESEARCH USE ONLY!**

SEXUALLY TRANSMITTED INFECTION PANEL

CT/NG/TV/MG Multiplex Positive Swab

Formulated with inactivated cells that contain whole genome *Chlamydia trachomatis* (CT), *Neisseria gonorrhoeae* (NG), *Trichomonas vaginalis* (TV), and human fibroblast cells

**MICROBIX BRAND
CATALOGUE NUMBER**

**For RUO: VP-S-12-M4
For IVD: RED-S-12-M4***

*Only available for commercial sale in USA and Canada

