

# Genital Ulcer Disease Panel

## QAPs™ for Molecular and Immunoassay Diagnostics

- Formulated with the target pathogen's whole genome
- Whole workflow controls (including extraction)
- Inactivated to meet clinical laboratory biosafety requirements
- Commutable and cross-platform compatible
- Liquid formulations are stable for 24 weeks when stored at 2-8°C; 1mL volume
- Copan FLOQSwab® formulations are stable for 1-2 years when stored at 2-30°C
- Verified performance on leading IVD platforms

### Compatible With

Product	Format	PTDx™	PROCEEDx™	REDx™	MDx	Immuno
<b>Single Analyte Products</b>						
HSV 1 Positive	Vial	PT-02-M1	VP-02-M1	Coming soon	✓	✓*
HSV 2 Positive	Vial	PT-23-M1	VP-23-M1	Coming soon	✓	✓*
HSV 1 Positive Swab	Swab	PT-S-02-01	Coming soon	Coming soon	✓	✓*
HSV 2 Positive Swab	Swab	PT-S-23-01	Coming soon	Coming soon	✓	✓*
<i>T.pallidum</i> Positive Swab	Swab	PT-S-80-01	Coming soon	Coming soon	✓	X
Varicella Zoster Virus Positive Swab	Swab	PT-S-03-01	Coming soon	Coming soon	✓	✓*
Monkeypox Congo Basin Positive Swab	Swab	Coming soon	Coming soon	Coming soon	✓	X
Monkeypox West African Positive Swab	Swab	Coming soon	Coming soon	Coming soon	✓	X
<b>Multiple Analyte Products</b>						
HSV1&2/VZV/Syphilis Positive Swab	Swab	PT-S-03-M1	VP-S-03-M1	RED-S-03-M1	✓	X
HSV 1 & 2 Positive	Vial	PT-02-M2	VP-02-M2	Coming soon	✓	✓*

Note: \*Products are QC- released with molecular assays but are compatible with immunoassays

## Intended Use

- REDx™ products are unassayed controls used for evaluating laboratory testing performance, procedures and workflows
- PROCEEDx™ products are used for assay verification/validation, training, and R&D purposes. Products are for RESEARCH USE ONLY!
- PTDx™ products are used as challenge samples in laboratory proficiency testing (PT) and external quality assessment (EQA) programs



## Need to Know More or Want to Order?

**Web:** [microbix.com](http://microbix.com) **Email:** [customer.service@microbix.com](mailto:customer.service@microbix.com)

**Phone:** 800-794-6694 toll free in North America or +1-905-361-8910 worldwide



**Microbix Biosystems Inc.**  
265 Watline Ave. Mississauga,  
ON Canada L4Z 1P3  
+1-905-361-8910  
[www.microbix.com](http://www.microbix.com)

[www.microbix.com](http://www.microbix.com)