

QAPs™ Negative Panel

QAPs™ for Molecular and Immunoassay Diagnostics



- Formulations do not contain infectious disease assay targets
- Whole workflow controls (including extraction)
- Inactivated to meet clinical laboratory biosafety requirements
- Commutable and cross-platform compatible
- Liquid formulations are stable for 1+ years when stored at 2-8°C; 0.5 mL and 1mL volumes
- Copan FLOQSwab® formulations are stable for 2 years when stored at 2-30°C
- Verified performance on leading IVD platforms

Product	Format	PTDx™	PROCEEDx™	REDx™	Human DNA Present	Compatible With	
						MDx	Immuno
STI Negative	Vial	PT-99-M1	VP-99-M1	RED-99-M1	Contains human DNA; compatible with most assays that detect specimen adequacy /cellularity with human housekeeping genes	✓	X
STI Negative	Vial	PT-99-M2	VP-99-M2	—	Contains human DNA; compatible with most assays that detect specimen adequacy /cellularity with human housekeeping genes	✓	✓
Negative	Vial	VP-99-01	VP-99-01	—	Does not contain human DNA	Only if MDx assay does not SAC requirements	✓
Respiratory Negative Swab (REDx™FLOQ®) Negative Swab (PTDx™ & PROCEEDx™FLOQ®)	Swab	PT-S-99-01	VP-S-99-01	RED-S-99-01	Does not contain human DNA	Only if MDx assay does not SAC requirements	✓
SARS-CoV-2 Negative	Vial	PT-99-M3	VP-99-M3	RED-99-M3	Contains human DNA in quantities that meet Sample Adequacy Control requirements for SARS-CoV-2 specific assays	✓	X
SARS-CoV-2 Negative Swab (REDx™FLOQ®) NAAT Negative Swab (PTDx™ & PROCEEDx™FLOQ®)	Swab	PT-S-99-M4	VP-S-99-M4	RED-S-99-M4		✓	X

Intended Use

- REDx™ products are unassayed controls that can be used for ongoing External Quality Control (QC), benchmarking laboratory testing performance, operator training & competency and end-to-end workflow evaluation.
- PROCEEDx™ products are used for assay verification/validation, training, and R&D purposes. Products are for RESEARCH USE ONLY!
- PTDx™ products are used as challenge samples in laboratory proficiency testing (PT) and external quality assessment (EQA) programs

Need to Know More or Want to Order?

Web: microbix.com **Email:** customer.service@microbix.com
Phone: 800-794-6694 toll free in North America or +1-905-361-8910 worldwide



Microbix Biosystems Inc.
 265 Watline Ave. Mississauga,
 ON Canada L4Z 1P3
 +1-905-361-8910
www.microbix.com