


Sexually Transmitted Infections

LIQUID FORMULATIONS



QAPs™ for STI Molecular and Immunoassay Diagnostics



- Contain target pathogen's whole genome
- Whole-process samples (including extraction)
- Inactivated to meet clinical laboratory biosafety requirements
- Commutable and cross-platform compatible
- Stable at 2-8°C; 1mL vials
- Verified performance on leading IVD platforms

Product	REDx™	PROCEEDx™	PTDx™	Compatible With	
				Molecular	Immuno
Single Analyte Products					
Chlamydia trachomatis Positive		VP-12-M1	PT-12-M1	✓	X
Chlamydia trachomatis Positive		VP-12-01	PT-12-01	✓	✓
Neisseria gonorrhoeae Positive		VP-17-M1	PT-17-M1	✓	X
Mycoplasma genitalium Positive	RED-63-01	VP-63-01	PT-63-01	✓	X
Trichomonas vaginalis Positive		VP-61-01	PT-61-01	✓	✓
Herpes Simplex Virus 1 Positive		VP-02-M1	PT-02-M1	✓	✓
Herpes Simplex Virus 2 Positive		VP-23-M1	PT-23-M1	✓	✓
Multiple Analyte Products					
Chlamydia trachomatis + Neisseria gonorrhoeae Positive		VP-12-M2	PT-12-M2	✓	X
Herpes Simplex Virus 1 & 2 Positive		VP-02-M2	PT-02-M2	✓	✓
Negative Analyte Products					
STI Negative	RED-99-M1	VP-99-M1	PT-99-M1	✓	X
STI Negative		VP-99-M2	PT-99-M2	✓	✓

Intended Use



- REDx™ products are unassayed controls used for evaluating laboratory testing performance, procedures and workflow
- PROCEEDx™ products are used for assay verification/validation, training, and R&D purposes. Products are for RESEARCH USE ONLY!
- PTDx™ products are used as challenge samples in laboratory proficiency testing (PT) and external quality assessment (EQA) programs

Need to Know More or Want to Order?

Web: microbix.com Email: customer.service@microbix.com

Phone: 800-794-6694 toll free in North America or +1-905-361-8910 worldwide

