


Viral Respiratory Panel

LIQUID FORMULATIONS



QAPs™ for Respiratory Molecular and Immunoassay Diagnostics



- Contain target pathogen's whole-viral genome
- Whole-process samples (including extraction)
- Inactivated to meet clinical laboratory biosafety requirements
- Commutable and cross-platform compatible
- Stable at **2-8°C**; 1mL vials
- Verified performance on leading IVD platforms

| Product | PROCEEDx™ | PTDx™ | Compatible With | |
|--|-----------|----------|-----------------|--------|
| | | | Molecular | Immuno |
| Single Analyte Products | | | | |
| Adenovirus Positive | VP-15-01 | PT-15-01 | ✓ | ✓ |
| Influenza A Positive | VP-13-01 | PT-13-01 | ✓ | ✓ |
| Influenza B Positive | VP-14-01 | PT-14-01 | ✓ | ✓ |
| Parainfluenza 3 Positive | VP-10-M1 | PT-10-M1 | ✓ | X |
| Respiratory Syncytial Virus (RSV) Positive | VP-07-01 | PT-07-01 | ✓ | ✓ |
| Multiple Analyte Products | | | | |
| Influenza A + Adenovirus Mixed Positive | VP-13-04 | PT-13-04 | ✓ | ✓ |
| Influenza B + Adenovirus Mixed Positive | VP-14-03 | PT-14-03 | ✓ | ✓ |
| Influenza A + Influenza B Mixed Positive | VP-13-02 | PT-13-02 | ✓ | ✓ |
| Influenza A + RSV Mixed Positive | VP-13-03 | PT-13-03 | ✓ | ✓ |
| Influenza B + RSV Mixed Positive | VP-14-02 | PT-14-02 | ✓ | ✓ |
| Adenovirus + RSV Mixed Positive | VP-07-02 | PT-07-02 | ✓ | ✓ |
| Negative Analyte Products | | | | |
| Respiratory Negative | VP-99-01 | PT-99-01 | ✓ | ✓ |

Refer to SARS-CoV-2 Product Brochures for an extensive list of SARS-CoV-2 Products

Intended Use

- PROCEEDx™ products are used for assay verification/validation, training, and R&D purposes. Products are for RESEARCH USE ONLY!
- PTDx™ products are used as challenge samples in laboratory proficiency testing (PT) and external quality assessment (EQA) programs

Need to Know More or Want to Order?

Web: microbix.com **Email:** customer.service@microbix.com

Phone: 800-794-6694 toll free in North America or +1-905-361-8910 worldwide

PROCEED PTD



Microbix Biosystems Inc.
265 Watline Ave, Mississauga,
ON Canada L4Z 1P3 +1-905-361-8910
www.microbix.com

www.microbix.com