

MICROBIX'S MOST POPULAR MULTIPLEX QUALITY CONTROLS

Enhance test monitoring and Performance with Microbix's
Human papillomavirus External Quality Controls
(QAPs™)



FEATURES

- **Whole workflow** controls
- **Cross-platform compatible**
- **Copan FLOQSwab®** stable for **2 years**
- Room-temperature stable swabs at **2-30°C**.
- Samples eluted in your choice of media
- **Verified performance** on IVD platforms
- Inactivated to **meet clinical laboratory biosafety requirements**

PRODUCT INTENDED USE

- **REDx™ products** are unassayed controls that can be used for ongoing **External Quality Control (QC)**, benchmarking **laboratory testing performance**, operator **training & competency** and end-to-end **workflow evaluation**.
- **PROCEEDx™ products** are used for **assay verification/validation, training, and R&D purposes**. Products are for **RESEARCH USE ONLY!**

HPV PANEL

HPV 16/18/45 Positive Swab

Formulated with inactivated human epithelial cells with integrated and episomal HPV 16, 18, and 45 genetic material desiccated on a Copan FLOQSwab®.

**MICROBIX BRAND
CATALOGUE NUMBER**

**For RUO: VP-S-62-M1
For IVD: Coming Soon**