

# MICROBIX'S MOST POPULAR MULTIPLEX QUALITY CONTROLS

Enhance test monitoring and Performance with Microbix's  
**STI External Quality Controls (QAPs™)**



## FEATURES

- **Whole workflow** controls
- **Cross-platform compatible**
- **Copan FLOQSwab®** stable for **2 years**
- Room-temperature stable swabs at **2-30°C**.
- Samples eluted in your choice of media
- **Verified performance** on IVD platforms
- Inactivated to **meet clinical laboratory biosafety requirements**

## PRODUCT INTENDED USE

- **REDx™ products** are unassayed controls that can be used for ongoing **External Quality Control (QC)**, benchmarking **laboratory testing performance**, operator **training & competency** and end-to-end **workflow evaluation**.
- **PROCEEDx™ products** are used for **assay verification/validation, training, and R&D purposes**. Products are for **RESEARCH USE ONLY!**

## SEXUALLY TRANSMITTED INFECTION PANEL

### CT/NG/TV/MG Multiplex Positive Swab

Formulated with inactivated Chlamydia trachomatis, Neisseria gonorrhoeae, Trichomonas vaginalis, Mycoplasma genitalium whole genome nucleic acid, and human cells, desiccated on a Copan FLOQSwab®.

**MICROBIX BRAND  
CATALOGUE NUMBER**

**For RUO: VP-S-12-M5**  
**For IVD: RED-S-12-M5\***

\*Only available for commercial sale in USA and Canada