


CONCLUSIONS

1. Full genome, quantified H5N1 cDNA control was developed by synthetic biology approach
2. Control tested in 75 labs, using 16 different qPCR assays
3. Can be used for typing assay development, assay verification/validation and EQA schemes

4. Control. Is non-infectious, stable and available in swab and liquid format
5. Microbix’s Platform approach can be utilized in pandemic and health emergencies as a rapid response to provide reference materials

AKNOWLEDGMENTS

We acknowledge the following EQA partners for providing the data used in this poster.



Pavel Zhelev, Director, Product Management
Email: pavel.zhelev@microbix.com



Download poster by scanning this QR code

