

Quality Assessment Products (QAPs)

www.microbix.com

PROCEEDx HPV Positive samples are designed to meet the full needs of the modern molecular diagnostic industry.

Advantages:

- Contain natural HPV types and human cells.
- 100% clinical sample equivalence.
- Prepared in HOLOGIC® ThinPrep® PreservCyt® STM.
- Inactivated to meet Clinical Lab Biosafety requirements.
- Designed to meet the instrument/device Sample Adequacy requirements.
- Industry leading panel representing a significant portion of hrHPV types.
- Dedicated hr (high-risk) Negative HPV 67 Positive sample for ruling out a test cross-reactivity.

Performance:

- PROCEEDx HPV Positive samples are compatible with Roche™ COBAS®, Cepheid® GeneXpert® and HOLOGIC® Aptima® platforms.
- Suitable for verification/validation on HPV DNA and RNA (TMA) testing systems.

Use:

- For Research Use Only (RUO).
- For instrument/device verification/validation and training.



**Need to know more
or want to order?**

visit us on the web at: microbix.com

email: customer.service@microbix.com

phone: 1-800-794-6694 toll free in North
America or +1-905-361-8910 worldwide



265 Watline Ave.
Mississauga, ON Canada L4Z 1P3
+1-905-361-8910
www.microbix.com

