



PRODUCT INFORMATION

Product Name:
PROCEEDx *Chlamydia trachomatis* Positive Sample

Catalogue Number:
VP-12-01

Product Description:
This product contains inactivated *Chlamydia trachomatis* in storage buffer.

SPECIFICATIONS

Strain:
Chlamydia trachomatis LGV type II strain 434

Storage:
Store this preparation at 2 °C to 8 °C.

Shelf Life:
72 weeks from manufacturing date.

Volume:
1 mL

USE

Instructions For Use:
Use PROCEEDx product as per sample handling instructions specified by device or instrument manufacturer.

Intended Use:
The product is to be used **for immunoassay** instrument (device) verification/validation, manufacturing QMS, educational and R&D purposes.

NOT FOR USE IN DIAGNOSTIC PROCEDURES. FOR RESEARCH USE ONLY.

ASSISTANCE

If you have any questions regarding the production, testing or use of this antigen, please send them together with any relevant data to Microbix Technical Services, at customer.service@microbix.com or fax 905-361-8911. Your complete satisfaction with the performance of this product is important to us.

HAZARDS




We are aware of no specific biological hazards associated with this product. The preparation has been inactivated and should contain no infectious material however it should be handled as if potentially infectious. Generally accepted good laboratory practices appropriate to biological reagents should be employed. The buffer solution contains 0.02% sodium azide (CAS No. 26628-22-8) as a preservative.

SPECIAL HANDLING

No special handling required for this product.

QUALITY CONTROL INFORMATION

| Test Kit | Target | Result |
|---|--------------------------------|--|
| Clear View Chlamydia MF Catalogue # 504524 | Chlamydia | Chlamydia Positive |
| QuickVue Chlamydia Catalogue # 0B006 | Chlamydia | Chlamydia Positive |
| Quidel Sofia RSV Catalogue # 20242 | Respiratory Syncytial Virus | RSV Negative |
| Quidel Sofia Influenza A & B Catalogue # 20254 | Influenza A & B | Influenza A Negative Influenza B Negative |
| SAS Adeno Catalogue # 067020 | Adenovirus | Adenovirus Negative |

| | |
|---|------------------------|
| REF | Catalog Number |
| LOT | Lot Number |
|  | Temperature Limitation |
|  | Expiration Date |
|  | Biological Risk |

