



## PRODUCT INFORMATION

**Product Name:**

PROCEEDx Adenovirus and Rotavirus Positive Sample

**Catalogue Number:**

VP-15-02

**Product Description:**

This product contains inactivated adenovirus and rotavirus in storage buffer.

## SPECIFICATIONS

**Strains:**

*Human adenovirus 2, strain Adenoid 6*, ICTV 00.001.0.01.010.00.102.006

*Simian rotavirus SA11*, ICTV 00.060.0.03.001.00.201.001

**Storage:**

Store this preparation at 2 °C to 8 °C.

**Shelf Life:**

72 weeks from manufacturing date.

**Volume:**

1 mL

## USE

**Instructions For Use:**

Use PROCEEDx product as per sample handling instructions specified by device or instrument manufacturer.

**Intended Use:**

The product is to be used **for immunoassay** instrument (device) verification/validation, manufacturing QMS, educational and R&D purposes.

**NOT FOR USE IN DIAGNOSTIC PROCEDURES. FOR RESEARCH USE ONLY.**

### ASSISTANCE

If you have any questions regarding the production, testing or use of this antigen, please send them together with any relevant data to Microbix Technical Services, at [customer.service@microbix.com](mailto:customer.service@microbix.com) or fax 905-361-8911. Your complete satisfaction with the performance of this product is important to us.

### HAZARDS




We are aware of no specific biological hazards associated with this product. The preparation has been inactivated and should contain no infectious material however it should be handled as if potentially infectious. Generally accepted good laboratory practices appropriate to biological reagents should be employed. The buffer solution contains 0.02% sodium azide (CAS No. 26628-22-8) as a preservative.

### SPECIAL HANDLING

No special handling required for this product.

# QUALITY CONTROL INFORMATION

Test Kit	Target	Result
Combi-Strip (Rota and Adeno) Catalogue # C-1004	Rotavirus & adenovirus	Positive
SAS Adeno Test Catalogue # 067020	Adenovirus	Positive
Immuno card STAT Rotavirus Catalogue # 750030	Rotavirus	Positive
Xpect Rotavirus Catalogue # R24655	Rotavirus	Positive
Binax Now RSV Catalogue # 430-022	Respiratory Syncytial Virus	Negative
Binax Now Flu A & B Catalogue # 416-000	Influenza A & B	Negative

<b>REF</b>	Catalog Number
<b>LOT</b>	Lot Number
	Temperature Limitation
	Expiration Date
	Biological Risk

